

BLOOD MIXER (ROLLER MIXER)

“ISO 9001:2008 & ISO 13485:2003 ”
CERTIFIED COMPANY
“CE MARKED PRODUCTS”



Model No : BM - 22

Roller Size : 4 Nos. X 340mm Lg.

Mains Input : 230v, 50Hz / 110v, 60Hz, AC

RPM : 30 RPM / 36 RPM

Mixing Time : 15 to 20 Minutes

Weight : 2.4 Kg. Net

Capacity : Vaccutainer 12 Nos.(of - 4ml)
Or bulb of 5ml x 24 Nos



LABLINE EQUIPMENTS PVT. LTD.
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,
BARODA 390010, GUJARAT, INDIA.
PH. NO. +91 265 3927448 / 49 / 50
FAX NO. +91 265 2632151
E-MAIL : haresh@lablineindia.com

VISIT US AT :- www.lablineindia.com